NSK

Cleaning Adaptor

CLA-20

OPERATION MANUAL

For Prophy Handpiece
For Perio Handpiece

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⚠ CAUTION

EN

- Read this Operation Manual before use to fully understand the product functions and file for future reference.
- Do not use this product for anything other than its intended purpose. Use this product according to the intended purpose and operation procedure described in the operation manual.
- Users are responsible for the operation control, maintenance and continual inspections of the Cleaning Adaptor.
- Do not allow any impact on to the Cleaning Adaptor. Do not drop the product. Deformation may cause the product to fail during use.
- Do not disassemble or alter the Cleaning Adaptor.
- Be sure to check the Cleaning Adaptor for cleanliness and damage before use.
- Only NSK "Perio Handpiece / Prophy Handpiece of Varios Combi Pro series" should be attached to the Cleaning Adaptor. Do not attach any other device.
- This product can only be used for manual cleaning of the handpiece. Do not place in washer-disinfector (auto-cleaning) or sterilizer. Reprocessing using washer-disinfector or sterilizer may cause damage or discoloration of the Cleaning Adaptor.
- This Cleaning Adaptor can be reused according to this operation manual.
- This product is consumable. Perform periodic maintenance check.
- Repairs cannot be made by the user.
- Do not use the following liquids to wipe, immerse, or clean the product: strong/super acid water, strong acid/alkaline chemicals, chlorine-containing solutions, and solvents such as benzine or thinner.
- To ensure proper function and safe use of the product, perform periodic maintenance checks in accordance with the operation manual.
- For reprocessing method of the Perio Handpiece/ Prophy Handpiece, refer to Varios Combi Pro series OPERATION MANUAL.
- No special training is required for this product.
- For cleaning and disinfection, use disinfectant ethanol or disinfectant isopropanol that does not contain any additives (hereinafter referred to as disinfectant alcohol).
- Use of disinfectants other than disinfectant alcohol may cause discoloration or cracking.
- Use appropriate disinfectant alcohol according to the regulations and guidelines of each country.
- Please report any serious incident that may occur in connection with this device and its accessories to the manufacturer as well as to your national authority.
- Hold the Perio/Prophy handpiece and this product firmly when spraying.

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
Cleaning System Figure (Manual Cleaning Only)

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EN Clean the surface of the Cleaning Adaptor using a lint-free cloth dampened with disinfectant periodically. Make sure that there is no deformation, damage or discoloration every three months. If any abnormalities are found, replace with a new one immediately.

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

EN In order to avoid any health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose of the product for you.

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
MODEL	REF			Remarks
CLA-20	ZA11900108	ø22 x H44mm	37g	Cleaning Adaptor for Perio Handpiece / Prophy Handpiece.

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NSK Europe GmbH EC REP
Elly-Beinhorn-Str.8, 65760 Eschborn, Germany

Visit our website

Specifications are subject to change without notice.

2025-02-13 CACF0221 002S