OPERATION MANUAL

NSK

CE 0197

Micromotor

NBX/NBX N

Please read this Operation Manual carefully before use, and file for future reference

MADE IN JAPAN

OM-E0561E

Powerful Partners®

Application
This micromotor is to be used for dental treatment only.

User
Dentists, therapists, hygienists or staff authorised by a dentist.

Prohibition
This micromotor is designed exclusively for general dental treatment. Do not use for surgical treatment.

◆ Classifications of equipment
- Not suitable for use in the presence of flammable anaesthetic mixture with air or with oxygen or nitrous oxide:
 – EQUIPMENT not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.
- Mode of operation.
 – Intermittent operation.

Classification	Degree of Danger or Damage and Seriousness
⚠CAUTION	Indicates an instruction where minor to medium injury or physical damage may occur.

⚠ CAUTION

- When operating the motor always consider the safety of the patient.
- This motor cannot operate without a control unit.
- Do not allow any impact on to the motor or drop it as this may cause a malfunction.
- Check for vibration, noise, and overheating as well as handpiece connection outside the patient's oral cavity before use. If any abnormalities are found, stop using the motor immediately and contact dealer.
- Do not exceed the motor speed recommended by the handpiece manufacturer.
- Do not lubricate the motor.
- Do not sterilize motor and motor cord.
- This motor cannot be cleaned and disinfected with a Thermo-Disinfector.
- If the motor has not been used for a long period, check for noise, rotation, and overheating before use.
- Air Requirements: dry, free from contamination and oil. Use a compressor with a dry air system. Install an air filter if necessary. Blow out the lines before installation.
- Do not use this motor under stress for long time as this will cause overheating.
- Do not pull the Motor excessively as this may loosen or disconnect the motor cord and cause malfunction.

- The user shall be responsible for operation and maintenance.
- The operator is responsible for correct operational control, maintenance and inspection.
- Surface temperature shall be more than 60°C without using a coolant air, to avoid this occurrence, be sure to use sufficient coolant air.
- When using the micromotor for more than 2 minutes continuously it will become hot. (The maximum surface temperature might reach 45℃) In such cases, stop using the micromotor until the motor has cooled down sufficiently.

1. Specification

Model		NBX	NBX N
Maximum Rotation Speed		40,000min^{-1}(rpm)	
Dimensions		Dø22 x L76 mm	
Input Voltage		DC 0 - 24V	
Spray Supplying water		More than 65 mL/min (Range: 0.2-0.25MPa)	
Chip Air		More than 1.5 L/min (Range: 0.2-0.25MPa)	
Coolant Air		More than 6.5 NL/min (Range: 0.2-0.25MPa)	
Intermittent Operation		3A for 2 min. , Air cooling for 5 min. (According to ISO60034-1 Type S3)	
Optics		White LED	–
LED	Voltage	DC 3.3V ± 0.5V	–
	Activate Voltage	2.5 - 4.0 V	–
	Consumption Current	0.2A Typ. (3.3V)	–
Use Environment		Temperature: 10 - 40˚C, Humidity: 30 - 75%RH*, Atmospheric pressure: 700 - 1,060hPa	
Transportation and Storage Environment		Temperature: -10 - 50˚C, Humidity: 10 – 85%RH*, Atmospheric pressure: 500 - 1,060hPa	

*No condensing allowed

2. Connecting/Disconnecting the Motor and the Motor Cord

Align and insert firmly the motor pins into the pin holes of the motor cord connector and fasten the motor cord nut securely. (Fig. 1)
To remove the motor cord from the micromotor unscrew and detach the motor cord nut and gently pull the micromotor and the motor cord connector apart.


Fig. 1

⚠ CAUTION

Make sure that the motor cord nut is firmly tightened. Loose connection will cause a water leak and motor may not function.

3. Connecting/Disconnecting the motor and the handpiece (Option)

To attach the E-type handpiece to the E-Type motor joint, align handpiece and motor, gently push handpiece on to the motor joint and then turn the handpiece until the LED light source in the motor line up in the correct position.

To remove the handpiece, gently pull the handpiece off the motor. (Fig. 2)

*Handpiece is not included.


Fig. 2

⚠ CAUTION

Care should be taken when using any handpiece attachment which requires lubrication. Following lubrication stand the handpiece vertically allowing the excess lubricant to drain and wipe the handpiece before attaching it to the motor.

4. Care

After each use, clean the surface of the motor with a soft cloth or nylon brush etc. (Do not use the metal brush), and wipe it with a moist soft cloth with surgical spirit.

⚠ CAUTION

- Do not autoclave (or any other high temperature sterilization) the motor.
- Do not lubricate the motor.
- Do not wash, soak or wipe the motor with/in potentially oxidizing solutions (strong acid, superacid solution) or sterilizing solution.
- Make sure the motor is kept in suitable atmospheric pressure, temperature, humidity and ventilation. Keep away from dust, salt and sulphur.

5. Replacing O-Rings

Replace the O-Rings when they look worn, as this may make it difficult to connect the handpiece or cause air and/or water leaks.
Remove the O-Rings from the E-Type joint with a blunt pointed tool and mount the new O-Rings into the groove. (Fig. 3)

*O-Ring (Blue) : Order Code 0313084070
O-Ring (Black) : Order Code 0312074080


Fig. 3

⚠ CAUTION

- There are 4 O-Rings on the E-Type motor joint. The blue O-Ring is thinner than the 3 black O-Rings. Make sure to replace all 4 O-Rings. (Fig. 4)
- Worn O-Rings may cause:
 - Air/Water leak
 - Air/Water blockages
 - Vibration
 - Difficulty connecting/disconnecting the handpiece
- Use ONLY original NSK O-Rings. If not, it may cause water and air leak.


Fig. 4

6. Troubleshooting

If you experience problems, please check the following again before consulting your dealer. If none of these is applicable and you are not able to resolve the problem a failure of this product is suspected and you should contact dealer.

Problem	Cause	Solution
Motor does not run	Motor Cord is not connected correctly	Check the connection
LED does not light	Reached the end of life expectancy	Contact dealer
Water leakage	Motor Cord is not connected correctly	Check the connection
	Water is leaking from the E-Type motor joint	Replace with new O-Ring

Symbols

Manufacturer

CE 0197 This conforms to CE European Directive of “Medical equipment directive 93/42/EEC.”

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CE 0197 The EU directive 93/42/EEC was applied in the design and production of this medical device.

NAKANISHI INC. 700 Shimohinata, Kanuma, Tochigi 322-8666, Japan www.nsk-inc.com

NSK Europe GmbH EC REP
Elly-Beinhorn-Strasse 8
65760 Eschborn
Germany

NSK France SAS
32 rue de Lisbonne
75008 Paris
France

NSK United Kingdom Ltd
Office 5, Gateway 1000
Arlington Business Park, Whittle Way
Stevenage, SG1 2FP, UK

NSK Dental Spain SA
C/ Módena, 43 El Soho-Európolis
28232 Las Rozas, Madrid
Spain

NSK America Corp
1800 Global Parkway
Hoffman Estates, IL 60192
USA

NSK America Latina Ltda
Dr. Placido Gomes, 610-Sala 104
A. Garibaldi
Joinville SC 89202-050, Brazil

NSK Oceania Pty Ltd
Unit 22, 198-222 Young St.
Waterloo, Sydney, NSW 2017
Australia

NSK Asia
1 Maritime Square
#09-33 HarbourFront Centre
099253, Singapore

NSK Middle East
Room 6EA701, 7th Floor, East Wing No.6
Dubai Airport Free Zone
PO Box 54316, Dubai, UAE

Specifications are subject to change without notice. 2015.01.20 01 Ⓜ

EMC INFORMATION

OM-E0569E 002

ENGLISH

Guidance and manufacturer' s declaration – electromagnetic emissions

The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.

Emissions test	Compliance test	Electromagnetic environment – guidance
RF emissions EN 55011 CISPR 11	Group 1	The product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions EN 55011 CISPR 11	Class B	The product is suitable for use in all establishments, including domestic establishments and those directly connected to the public low–voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions EN/IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions EN/IEC 61000-3-3	Complies	

EMC Measurement Conditions

No.	Interface(s)	Max. cable length, shielding	Cable classifications
1	AC input (with plug)	Un–Specified	AC Power Supply Line
2	AC output	5.0m, Shielded	AC Power Supply Line
3	Motor line	2.2m, Un–shielded	DC Power Line
4	Lamp line	2.3m, Un–shielded	DC Power Line

Guidance and manufacturer' s declaration – electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) EN /IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst EN/IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ±1kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN/IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line to line ± 2 kV lines to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN/IEC 61000-4-11	< 5 % Ut (> 95 % dip in Ut) for 0.5 cycle 40 % Ut (60 % dip in Ut) for 5 cycles 70 % Ut (30 % dip in Ut) for 25 cycles < 5 % Ut (> 95 % dip in Ut) for 5 sec	< 5 % Ut (> 95 % dip in Ut) for 0.5 cycle 40 % Ut (60 % dip in Ut) for 5 cycles 70 % Ut (30 % dip in Ut) for 25 cycles < 5 % Ut (> 95 % dip in Ut) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the product requires continued operation during power mains interruptions, it is recommended that the product be powered from an uninterruptible power supply or a battery.
Power frequency (50 / 60 Hz) magnetic field EN/IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE : Ut is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer' s declaration – electromagnetic immunity

The product is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the product, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance. Where $\sqrt{P}$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey[a], should be less than the compliance level in each frequency range[b]. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF EN/IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	$d = 1.2\sqrt{P}$
Radiated RF EN/IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz

NOTE 1 At 80MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobiles radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters an electromagnetic site survey should be considered. If the measured field strength in the location in which the product is used exceeds the applicable RF compliance level above, the product should be observed to verify normal operation. If abnormal performance is observed additional measures may be necessary, such as reorienting or relocating the product.

b Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the product

The product is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the product as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of the transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz $d = 1.2\sqrt{P}$	80MHz to 800MHz $d = 1.2\sqrt{P}$	800MHz to 2.5GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations, Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

 Marking on the outside of Equipment or Equipment parts that include RF transmitters or that apply RF electromagnetic energy for diagnosis or treatment.

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NSK Asia Pte Ltd