NSK

VarioSurg 3 Handpiece

REPROCESSING MANUAL

VS3-LED-HPSC

CE 0197 MADE IN JAPAN

EN This reprocessing manual describes the procedure of reprocessing the VarioSurg 3 Handpiece. Carefully read the operation manual of the Ultrasonic System for other procedures.

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General Information

EN User: Qualified professionals (dentist, oral surgeon)

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EN Intended purpose: Dentistry and oral surgery (implant site preparation, bone cutting, maxillary sinus floor (mucus membrane) elevation, osteoplasty, bone resection in radectomy, periodontal operation, prosthesis maintenance, surgical endodontic procedure.

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EN Immediately after treatment, the handpiece should be cleaned, disinfect, sterilized and stored in accordance with this reprocessing manual.

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WARNING

EN • Failure to properly maintain the handpiece as per below may cause infection, product failure or overheating leading to burn injuries. • Be sure to always carefully read the operation manual for the Ultrasonic System before use to understand the usage and performance of the instrument.

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CAUTION

EN • Use only cleaning disinfectant products listed in this manual. • Follow any additional local directives, standards and guidelines for cleaning, disinfection, and sterilization. • Follow local rules, regulations and guidelines regarding the reprocessing of devices. • Handpiece is delivered in a non-sterile condition and requires reprocessing prior to first use and after each patient. • Be sure to clean the handpiece immediately (within 30 minutes) after use to remove residue. • The Handpiece can not be removed from the Handpiece Cord. • Do not clean, immerse, or wipe with electrolyzed-oxidizing water (strong acidic water, super acidic water), strong acidic and strong alkaline liquid medicine, solvent containing chlorine, benzene, or thinner. • Use clean and dry compressed air for cleaning. • When using in a washer-disinfector, ensure that the Handpiece and Handpiece Cord are stable by placing them in a basket or container. • Handpiece must be removed from the washer-disinfector immediately (within 1 hour) after the cleaning, disinfecting and drying cycle is complete to prevent corrosion. • After using washer- disinfector, dry the product completely. The remaining moisture causes internal corrosion etc. • Do not autoclave the handpiece with other instruments even when it is in a pouch. • This is to prevent possible discoloration and damage to the product from chemical residue on other instruments. • Do not heat or cool the product too quickly. Rapid change in temperature could cause damage to the product. • Do not set the sterilization temperature parameter higher than 135°C. • Contact the sterilizer manufacturer for detailed information about cycle temperatures. • Clean the handpiece prior to disinfection. If blood remains on the internal surface it can become clotted and cause product failure. • For cleaning and disinfection, use disinfectant ethanol or disinfectant isopropanol that does not contain any additives (hereinafter referred to as disinfectant alcohol). Use of disinfectants other than disinfectant alcohol may cause discoloration or cracking. • Immediately after sterilization is complete (within 1 hour), remove the handpiece from the

sterilizer. Failure to do so may cause corrosion. • To ensure proper function and safe use of the medical device, perform periodical maintenance checks in accordance with the operation manual of the ultrasonic system. • No special training is required for this device. • For sterilization method, refer to this operation manual. • Do not clean the handpiece in an ultrasonic cleaner. • Handpiece is a consumable item. • Repeated autoclaving may cause the handpiece to become discolored due to heat. However, this is due to properties of the product and is not problem in terms of quality. • When the product has been reprocessed 250 times or has been in use for one year, it is recommended that a periodic inspection be performed by an authorized NSK service. • The product lifetime and product' s ability to operate correctly are mainly determined by mechanical stresses during use and the chemical influences due to processing. If you detect any abnormality such as heat generation, vibration or noise, or if the exterior markings (serial number or model number) have become unreadable before or during use, stop using the product immediately and request your nearest NSK dealer, or one of the NSK branches for a periodic inspection. • Report all serious events related to the product immediately, to the manufacturer and responsible national authority.

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

EN Wait for activation time following the instructions provided by the cleaning/disinfectant manufacturer.

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
EN Use appropriate disinfectant alcohol according to local regulation and guidelines in force.

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
EN In order to avoid any health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose of the product for you.

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
EN * Avoid any risk of freezing of irrigation solution When operated outside the above range, there is a risk of failure.

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WARNING

EN WARNING: Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.

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CAUTION

EN CAUTION: Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.

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NSK Europe GmbH EC REP