NSK

OM-C0591E 003

OPERATION MANUAL

CE 0197

MADE IN JAPAN

Ti-Max Z　Ti-Max X
Ti-Max Z45L　S-Max M
Ti-Max nano

Download operation manuals

www.nsk-dental.com/support/userguide.html

NAKANISHI INC.　WWW.nsk-dental.com
700 Shimohinata, Kanuma, Tochigi 322-8666, Japan

NSK Europe GmbH　EC REP
Elly-Beinhorn-Str. 8, 65760 Eschborn, Germany

EN

- User: Qualified professionals (dentist)　Intended purpose: Z15L / Z24L / Z25L / Z85L / Z95L / Z99L / nano15LS / nano25LS / nano95LS / nano65LS / X12L / X12 / X15L / X15 / X20L / X25L / X25 / X85L / X85 / X95L / X95 / X95EX / X65L / X65 / M15L / M15 / M25L / M25 / M94L / M95L / M95 / M65: The dental handpiece is intended for the following application(s): Caries removal, cavity and crown preparation, removal of dental restorations (fillings and prosthesis), finishing and polishing of teeth and dental restorations (preparation and adjustment), and removal of third molars. Patient population: Any patient (ages 2 and over). The user's discretion is required prior to use.

- Depressing the Push Button while the handpiece is in rotation may lead to overheating, causing burn injuries or product failure. Avoid the push button to contact with any oral tissue. • Supply coolant water through the spray nozzle while using the handpiece. Other water supply methods than recommended may lead to overheating causing burn injuries, or product failure. • Do not use the handpiece for polishing purposes. Polishing paste could enter the handpiece, causing the failure of the Push Button or cause burn injuries due to the overheating of the handpiece head. • Keep any debris or other foreign materials away from inside the gear or the handpiece. Foreign materials remaining inside may lead to overheating, causing burn injuries or other accidents. • Ball bearings, in time, will become worn and will require replacement. Be sure to inspect the handpiece before use. If abnormal vibration or noise are found, the bearings or other internal parts may require replacement. Damaged parts may lead to handpiece overheating causing burn injuries. If overheating is suspected, stop using the handpiece immediately and contact your Authorized NSK Dealer.

- When operating the product always consider the safety of the patient. • Users are responsible for the operational control, maintenance and continual inspection of this product. • Do not disassemble or alter the handpiece and cartridge, etc. except as recommended by NSK in this Operation Manual. • Should the handpiece function abnormally during use, stop using the handpiece immediately and contact your Authorized NSK Dealer. • The product is delivered in a non-sterile condition. Be sure to sterilize the product prior to the first use or prior to use. • Follow the specified procedures in this operation manual as to the burr to be used, etc. to prevent accidental release or breakage of the burr during rotation. • Perform regular function and maintenance checks. • To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during surgery. • The operation of the handpiece is permitted only on dental units which correspond to the standards IEC 60601-1 (EN 60601-1) and IEC 60601-1-2 (EN 60601-1-2). • Contact your Authorized NSK Dealer if repairs are necessary. Z95L / Z96L / X95L / X95EX / M95L / M94L / M95 / nano95LS: • Use a surgical burr with its maximum length of 25mm or longer. Using such bur with this handpiece may lead to premature abrasion of the bearings. In addition, sudden release, bend or breakage of the bur may occur. Z94L / Z85L / X85 / X85EX / M85 / nano85LS: • Use a burr shank Ø2 bur for this handpiece. Using a bur with a shank Ø1 bur for this handpiece has the miniature head. Using a long bur with this handpiece may lead to premature abrasion of the bearings. In addition, sudden release, bend or breakage of the bur may occur. Z99L / Z45L: • When using a bur with a length of 20-25mm, the working diameter must be within the value specified as below. Diamond Burr: Ø2 mm or less, Carbide Bur: Ø1mm or less Be careful of overheating when using some types of carbide bur, as the coolant water (= "water jet" mode) cannot reach the bur tip because of the shape of the tip. X65L / X65 / M65 / nano65LS: • If the chuck is left closed for an extended period without a bur inserted, it could become difficult to insert the bur. • Be sure to insert the Bur Stopper in the chuck before inserting a CA bur. Once trapped in the chuck without the Bur Stopper, the bur will be difficult to remove. • Do not use a HP bur in combination with the Bur Stopper as the bur will not be fully secured and may disengage from the chuck, causing an injury. • Do not use Short Shank CA burs as the bur will not be fully secured and may disengage from the chuck, causing injury. • Please report any serious incident that may occur in connection with this product and its accessories to the manufacturer as well as to your national authority. • Do not exceed the Max. Rotation Speed (Motor) shown in the specifications in this operation manual. • Do not apply excess pressure to the bur as it may break or bend and become difficult to remove. • Always follow the instructions provided for the bur manufacturer regarding the bur usage and rotation speed, etc.

- Substances which are carcinogenic, mutagenic or toxic to reproduction (「CMR」), of category 1A or 1B, in accordance with Part 3 of Annex VI to Regulation (EC) No 1272/2008 of the European Parliament and of the Council, or substances having endocrine-disrupting properties for which there is scientific evidence of probable serious effects to human health and which are identified either in accordance with the procedure set out in Article 59 of Regulation (EC) No 1907/2006 of the European Parliament and of the Council or, once a delegated act has been adopted by the Commission pursuant to Article 57 of Regulation (EU) No 528/2012 of the European Parliament and the Council, in accordance with the criteria that are relevant to human health amongst the criteria established therein: lead (Pb), cobalt (Co). This product can be safely used in the treatment of children, pregnant or breastfeeding women, or other patient groups considered particularly vulnerable to such substances.

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