NSK

OPERATION MANUAL

Please read this OPERATION MANUAL carefully before use, and file for future reference.

E-Type Handpiece for Laboratory Work

VR-EB

CE MADE IN JAPAN

OM-H0228E 001

1. User and Indications for Use

Users: Dental Technicians and Dentists
Indications for Use: For use in Dental Laboratories for cutting and polishing

2. Precautions for handling and operation

- Please read these precautions carefully and use only as intended or instructed.
- Safety instructions are intended to avoid potential hazards that could result in personal injury or damage to the device. Safety instructions are classified as follows in accordance with the seriousness of the risk.

Class	Degree of Risk
⚠ WARNING	Hazard that could result in serious injury or damage to the device if the safety instructions are not correctly followed.
⚠ CAUTION	Hazard that could result in light or moderate injury or damage to the device if the safety instructions are not correctly followed.
NOTICE	General product specification information highlighted to avoid product malfunction and performance reduction.

⚠ WARNING

- This product is only for use in dentistry. Do not use it for medical treatment.
- Conform to the allowed rotation speed specified by the manufacturer/distributor of the bur. Failure to Do so may lead to injury due to the bur breaking and scattering or the Motor Handpiece moving violently.
- Do not use burs with instability, curving, damage, or wear on the shank area. Doing so may lead to injury due to the bur breaking or flying out, etc.
- When using a disk, use it at as low a rotation speed as possible to prevent it from scattering due to breaking, etc.
- Since a grindstone may be unbalanced even if it is new, make sure to perform dressing before using it. Using an unbalanced grindstone may lead to injury due to it breaking and scattering.

⚠ CAUTION

- The product must be used by a person who has received professional instruction as a health care provider, such as a dental technician or dental physician, in a dental clinic, hospital or other medical institution.
- Users are responsible for the operational control, maintenance and continual inspection of this product.

⚠ CAUTION

- Rotate the bur to check for vibrations, noise, or temperature (heat) being emitted. Should the Handpiece function abnormally at this time, stop using it immediately and contact your Authorized NSK Dealer. (Refer to “6. Checking the Handpiece Before Each Use”)
- Do not disassemble or alter the handpiece except as recommended by NSK in this Operation Manual.
- Do not allow any impact on to the product. Do not drop the product.
- Operators and all others in the area must wear eye protection, a mask, and gloves when operating the Handpiece.
- Operate the Handpiece at speeds under the maximum rotation speed of 35,000min^{-1}.
- When using a bur with a large cutting blade (4mm or larger in diameter), use it at as low a rotation speed as possible.
- Do not use burs with problems listed below as the bur may break, seize up or disengage from the chuck.
 – Bent, deformed, worn, rusted, broken, deficient bur.
 – Bur which is cracked on the edge or axis.
 – Non-ISO (EN ISO) standard, or tampered bur.
- Always keep the bur shank clean. Debris in the chuck could cause poor bur concentricity or low chuck retention force.
- Always insert the bur all the way into the chuck. If insertion is insufficient, premature failure of the handpiece head bearings, or accidental release of the bur, may occur.

⚠ CAUTION

- Always follow the instructions provided by the bur manufacturer.
- When a bur is rotating, do not turn the bur lock ring in the direction. Doing so will cause the Handpiece to break. Particularly when replacing a bur, do so after rotation has completely stopped.
- Should the handpiece function abnormally during use, stop using the handpiece immediately and contact your Authorized NSK Dealer.
- Clean the Handpiece on a timely basis. Follow the procedure in “7. Post-use Maintenance” in the Operation Manual to perform cleaning.
- Do not use the following fluids to wipe, immerse or clean the product; strong / super acid water, strong acid / alkaline chemicals, chlorine-containing solutions, solvents such as benzine or thinner.
- Do not lubricate the Motor Handpiece with oil. Doing so might cause it to generate heat or malfunction.
- Perform regular function and maintenance checks. (Refer to “8. Periodical Maintenance Checks”)
- If the handpiece has not been used for a long period, rotate the handpiece and check for noise, vibration and overheating before use.
- To avoid clinical downtime it is recommended that a spare be kept on hand in case of a breakdown during treatment.

NOTICE

- Keep a cutting bur or test bur attached even when the product is not in use.
- Contact your Authorized NSK Dealer if the product fails.

3. Accessory List

No.	Name	Quantity
1	Wrench for Adjusting Chuck	1

4. Connecting and Disconnecting the Handpiece from the Motor

4-1 Connecting

1) Insert the handpiece direct into the motor (Fig. 1). When connecting an optic handpiece, turn the handpiece until it clicks and locks into position.
2) Confirm that the handpiece is securely connected to the motor.

4-2 Disconnecting

Hold the motor and the handpiece, then pull apart.


Fig. 1

⚠ CAUTION

- Do not connect or disconnect the handpiece until the motor has completely stopped.
- Connect ONLY to E type motors (ISO 3964 (EN ISO 3964)).

5. Inserting and Removing the Bur

5-1 Removing the Bur

1) Turn the bur lock ring in the direction until it clicks to loosen the chuck.
2) Pull out the bur.

5-2 Inserting the Bur

1) Insert the bur in all the way.
2) Turn the bur lock ring in the direction until it clicks to tighten the chuck. After attaching the bur, pull it back and forth to confirm that it is connected securely.


Fig. 2

⚠ WARNING

- Conform to the allowed rotation speed specified by the manufacturer/distributor of the bur. Failure to Do so may lead to injury due to the bur breaking and scattering or the Motor Handpiece moving violently.
- Do not use burs with instability, curving, damage, or wear on the shank area. Doing so may lead to injury due to the bur breaking or flying out, etc.
- When using a disk, use it at as low a rotation speed as possible to prevent it from scattering due to breaking, etc.
- Since a grindstone may be unbalanced even if it is new, make sure to perform dressing before using it. Using an unbalanced grindstone may lead to injury due to it breaking and scattering.

⚠ CAUTION

- When using a bur with a large cutting blade (4 mm or larger in diameter), use it at as low a rotation speed as possible.
- Always keep the bur shank clean. Debris in the chuck could cause poor bur concentricity or low chuck retention force.
- Always insert the bur all the way into the chuck. If insertion is insufficient, premature failure of the handpiece head bearings, or accidental release of the bur, may occur.
- Always follow the instructions provided by the bur manufacturer.
- When a bur is rotating, do not turn the bur lock ring in the OPEN direction. Doing so will cause the Handpiece to break. Particularly when replacing a bur, do so after rotation has completely stopped.

6. Checking the Handpiece Before Each Use

Before using the product, rotate the bur to perform an inspection. If you notice a problem such as bur instability, vibrations, generation of heat, or abnormal sounds during an inspection or use of the product, cease operation immediately and contact your Authorized NSK Dealer.

7. Post-use Maintenance

Perform the following maintenance after use.

7-1 Preparation

1) Remove the bur.
2) Disconnect the handpiece from the motor.
3) Wipe the surface of the Handpiece with a damp cloth, and then wipe it with a cloth dampened with rubbing alcohol, etc.

7-2 Cleaning and Replacing the Chuck

7-2-1 Removing the Chuck

1) Turn the bur lock ring to the OPEN position.
2) Turn the chuck counterclockwise with the included wrench for adjusting chuck (Fig. 3).


Fig. 3

3) Remove the chuck (Fig. 4).

7-2-2 Cleaning the Chuck

In order to maintain the precision of the chuck, remove it regularly for cleaning with an ultrasonic cleaner, etc. It is recommended that you clean it at least once a week to ensure the product can be used safely.


Fig. 4

⚠ CAUTION

- Neglecting to clean the chuck for a long time is very dangerous because wax, gypsum, etc., will accumulate in the chuck, causing the bur to be insecure and causing runout.

7-2-3 Attaching the Chuck

1) Turn the bur lock ring to the OPEN position.
2) Push the chuck in all the way with a bur inserted.
3) Turn the chuck clockwise with your finger until it stops (Fig. 5).
 Then lock the bur lock ring to hold the bur securely.


Fig. 5

⚠ CAUTION

- After attaching the chuck, pull it back and forth to confirm that it is connected securely.

NOTICE

- Apply a thin layer of oil to the chuck before insertion.

8. Periodical Maintenance Checks

Every 3 months perform periodical maintenance checks, referring to the check sheet below. If any abnormalities are found, contact your Authorized NSK Dealer.

Points to check	Details
Rotation	Rotate the handpiece and check for abnormalities such as abnormal rotation, vibration, noise, and overheating.

9. Problems and Solutions

Problem	Cause	Solution
The bur does not rotate with the bur lock ring closed.	The ball bearings have burned out or have foreign material inside them.	Contact your Authorized NSK Dealer.*
Heat is emitted during rotation.	There is wear on the ball bearings due to foreign material inside them.	Contact your Authorized NSK Dealer.*
Vibrations/noise is emitted during rotation.	There is wear on the ball bearings due to foreign material inside them.	Contact your Authorized NSK Dealer.*
	A curved bur was used.	Replace the bur.
The bur has poor concentricity.	Foreign material may be attached to the inside of the chuck.	Clean the inside of the chuck.
	The chuck is worn.	Replace the chuck.
	The ball bearings are worn.	Contact your Authorized NSK Dealer.*

Problem	Cause	Solution
The bur disconnects.	The chuck is loose.	Tighten the chuck.

*The product cannot be repaired at a dental clinic, etc.

10. Specifications

Model	VR-EB
Max. Rotation Speed (Motor)	35,000min^{-1}
Max. Rotation Speed (Handpiece)	35,000min^{-1}
Bur Type	ISO 1797-1 (EN ISO 1797-1) Type2 ø2.35mm HP Bur
Use Environment	Temperature: 10 - 40°C, Humidity: 30 - 75% (No Condensation)
Transportation and Store Environment	Temperature -10 - 50°C, Humidity: 10 - 85%, Pressure: 500 - 1,060hPa

11. Symbol

Conforms to CE European Directive.

Manufacturer.

Authorized representative in the European community.

12. Warranty

NSK products are warranted against manufacturing errors and defects in materials. NSK reserves the right to analyze and determine the cause of any problem. Warranty is voided should the product be not used correctly or for the intended purpose or has been tampered with by unqualified personnel or has had non NSK parts installed. Replacement parts are available for seven years beyond discontinuation of the model.

13. Spare Parts List

Model	Order Code
Chuck-2.35A	H203A180A
Wrench for Adjusting Chuck	K015751

14. Disposing product

In order to avoid the health risks of operators handling the disposal of medical equipment, as well as the risks of environmental contamination caused thereof, a surgeon or a dentist is required to confirm the equipment is sterile. Ask specialist firms who are licensed to dispose of specially controlled industrial wastes, to dispose the product for you.

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NSK Europe GmbH EC REP